

Diameters and lengths - straight stents

GSS™ TD / GSS™ TF / GSS™ BD / DUMON® BB

Determine the desired stent length (stenosis length + 10 mm) and diameter and select from the available stents.



Custom-made stents are possible. Contact us!

Length (mm) ▶	20	30	40	50	60	70	80	90	100	110
Outer Diameter (mm) ▼	20	30	40	50	60	70	80	90	100	110
5	025301S20 (BB)	025301S30 (BB)	025301S40 (BB)	025301S50 (BB)						
6	026201S20 (BB)	026201S30 (BB)	026201S40 (BB)	026201S50 (BB)						
7	026501S20 (BB)	026501S30 (BB)	026501S40 (BB)	026501S50 (BB)						
8	026701S20 (BB)	026701S30 (BB)	026701S40 (BB)	026701S50 (BB)						
10	01BD1020	01BD1030	01BD1040	01BD1050	01BD1060	01BD1070				
11	01BD1120 01TD1120	01BD1130 01TD1130	01BD1140 01TD1140	01BD1150 01TD1150	01BD1160 01TD1160	01BD1170 01TD1170	01TD1180			
12	01BD1220 01TD1220	01BD1230 01TD1230 01TF1230	01BD1240 01TD1240 01TF1240	01BD1250 01TD1250 01TF1250	01BD1260 01TD1260 01TF1260	01BD1270 01TD1270 01TF1270	01BD1280 01TD1280 01TF1280			
13		01TD1330 01TF1330	01TD1340 01TF1340	01TD1350 01TF1350	01TD1360 01TF1360	01TD1370 01TF1370	01TD1380			
14		01TD1430 01TF1430	01TD1440 01TF1440	01TD1450 01TF1450	01TD1460 01TF1460	01TD1470 01TF1470	01TD1480			
15		01TD1530 01TF1530	01TD1540 01TF1540	01TD1550 01TF1550	01TD1560 01TF1560	01TD1570 01TF1570	01TD1580 01TF1580	01TD1590 01TF1590	01TD15100 01TF15100	01TD15110 01TF15110
16		01TD1630	01TD1640 01TF1640	01TD1650 01TF1650	01TD1660 01TF1660	01TD1670 01TF1670	01TD1680 01TF1680	01TD1690 01TF1690	01TD16100 01TF16100	01TD16110
18			01TD1840 01TF1840	01TD1850 01TF1850	01TF1860	01TF1870	01TF1880	01TF1890	01TF18100	01TF18110
20			01TF2040	01TF2050	01TF2060	01TF2070	01TF2080	01TF2090	01TF20100	01TF20110



GSS™ TD



GSS™ TF



GSS™ BD



DUMON® BB

